Traxtion® hexagonal snare

Anse hexagonale Traxtion®
Traxtion® hexagonale Schlinge
Traxtion® - ansa esagonale
Traxtion® - lazo hexagonal
Ansa hexagonal Traxtion®
Traxtion® heksagonal slyngje
Traxtion® zeshoekige poliepsnoerder
Traxtion® altigen snare
Traxtion® 6각형 올가미
Петля полипэктомическая
невращающаяся шестиугольная
Traxtion®

Reorder No. 00711113

Référence de commande
Nachbestell-Nr. 00711113
N. di riordino 00711113
N.º de pedido 00711113
N.º de encomenda 00711113
Genbestillingsnr. 00711113
Nabestelnr. 00711113
Yeni Sipariş Numarası: 00711113
주문 번호 00711113
№ для заказа 00711113

Short Throw™ snare – mini-oval

Anse Short Throw™ - mini-ovale
Short Throw™ Schlinge – Mini, oval
Short Throw™ - ansa – mini - ovale
Short Throw™ - lazo – mini-ovalado
Ansa Short Throw™ – mini-oval
Short Throw™ slyngje – mini-oval
Short Throw™ poliepsnoerder – mini-ovaal
Short Throw™ snare - mini-oval
Short Throw™ 올가미 - 소형 타원
Петля полипэктомическая
невращающаяся овальная размера
мини Short Throw™

Reorder No. 00711018

Référence de commande 00711018
Nachbestellnr. 00711018
N. di riordino 00711018
N.º de pedido 00711018
N.º de encomenda 00711018
Genbestillingsnr. 00711018
Nabestelnr. 00711018
Yeni Sipariş Numarası: 00711018
주문 번호 00711018
№ для заказа 00711018

Short Throw™ - snare – enteroscope

Short Throw™ - anse - entéroscope
Short Throw™ Schlinge – Enteroskop
Short Throw™ - ansa – enteroscopio
Short Throw™ - lazo – enteroscopia
Ansa Short Throw™ – enteroscópio
Short Throw™ slyngje – enteroskop
Short Throw™ - poliepsnoerder – enteroscoop
Short Throw™ - snare – enteroskop
Short Throw™ - 올가미 – 장경
Петля полипэктомическая
невращающаяся для использования
с энтероскопом

Reorder No. 00711101

Référence de commande 00711101
Nachbestell-Nr. 00711101
N. di riordino 00711101
N.º de pedido 00711101
N.º de encomenda 00711101
Genbestillingsnr. 00711101
Nabestelnr. 00711101
Yeni Sipariş Numarası: 00711101
주문 번호 00711101
№ для заказа 00711101

Short Throw™ - snare – oval

Short Throw™ - anse – ovale
Short Throw™ Schlinge – oval
Short Throw™ - ansa – ovale
Short Throw™ - lazo – ovalado
Ansa Short Throw™ – oval
Short Throw™ slyngje – oval
Short Throw™ - poliepsnoerder – ovaal
Short Throw™ - snare – oval
Short Throw™ - 올가미 – 타원
Петля полипэктомическая
невращающаяся овальная Short
Throw™

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2797

This product is not made with natural rubber latex.

Product	product number	sheath diameter (mm)	length (cm)	approximate snare size (cm)	active cord style
Traxtion® hexagonal snare	BX00711113	2.4	230	2.5 x 4.0	Olympus
Short Throw™ snare - mini-oval	BX00711018	2.4	230	1.25 x 3.0	Olympus
Short Throw™ snare - oval	BX00711019	2.4	230	2.5 x 5.4	Olympus
Short Throw™ snare - enteroscope	BX00711101	2.4	350	2.5 x 5.4	Olympus
Rotator™ snare – oval	BX00711120	2.4	230	2.5 x 5.4	Olympus
Rotator™ snare - mini-oval	BX00711121	2.4	230	1.25 x 3.0	Olympus

Intended Use:

The Rotator™ snares, Short Throw™ snares and Traxtion® hexagonal snare are electrosurgical devices designed to be used endoscopically to grasp, dissect and transect tissue during gastrointestinal (GI) endoscopic procedures. These devices are used with diathermic energy.

Warnings and Precautions:

- Consult the medical literature relative to techniques, complications and hazards prior to the performance of any endoscopic procedure.
- Endoscopic procedures should only be performed by persons having adequate training and familiarity with endoscopic techniques in snare polypectomy and injection with an endoscopic injection needle.
- A thorough understanding of the technical principles, clinical applications, and risks associated with monopolar diathermic energy and cold snaring (non-electrical) polypectomy and tissue resection is necessary before using this product.
- Endoscopic polypectomy with monopolar diathermic snares should not be performed without a thorough understanding of the principles of diathermic energy.
- Do not attempt to reuse, reprocess, refurbish, remanufacture, or resterilize this device. STERIS Endoscopy did not design this device nor is it intended to be reused, reprocessed, refurbished, remanufactured, or resterilized. Performing such activities on this disposable medical device presents a safety risk to patients (i.e. compromised device integrity, cross-contamination, infection).
- These devices are compatible with an endoscope accessory channel of 2.8mm or larger.
- Short strokes, 1"-1.5" (2.5cm – 3.8cm) in length, are recommended throughout device passage to avoid sheath kinking.
- This device is not designed for use in side viewing endoscopes.
- The following conditions may not allow a snare device to function properly:
 - (1) advancing the handle to the open position with too much speed or force,
 - (2) attempting to pass or open the device in an extremely articulated endoscope,
 - (3) attempting to actuate the device in an extremely coiled position and/or
 - (4) actuating the device when the handle is at an acute angle in relation to the sheath.
- If there is lack of visible cautery effect, consider the following:
 - (1) Check the active cord for secure connection to the snare handle and the electrosurgical generator.
 - (2) Ensure that the return electrode is secure and properly connected to the patient and the electrosurgical generator.
 - (3) Consult the electrosurgical generator manufacturer's instructions for use for proper settings and use of the generator.
 - (4) **In order to ensure that the insulating properties of the monopolar diathermic snares are not compromised, do not exceed the maximum rated peak voltage of 2500 for cut mode and 2500 for coagulation mode**

Caution - Output power setting of the electrosurgical unit should be set as low as possible to achieve the intended purpose. If the power setting of the generator is not known, set the power level at the minimum value of the expected range and cautiously increase power to achieve the desired tissue effect

- Observe for whitening of the tissue as visible evidence of electrocoagulation.
- User should apply caution during diathermic procedures to avoid cutting tissue before the application of cautery.
- To avoid patient injury, do not apply excessive force to the snare while resecting tissue.
- Disengage the snare from the polyp if there is a risk of complication.

Contraindications:

- Contraindications include, but are not limited to, those specific to any endoscopic procedure.
- Coagulopathy
- Poor bowel prep
- Uncooperative patient
- Consult the appropriate medical personnel regarding use of this device in patients with an implanted electronic device.

Prior to Use:

1. Prior to clinical use, inspect and familiarize yourself with the device. If there is evidence of damage, **do not use this product and contact your local Product Specialist.**
2. Ensure that the device is compatible with an endoscope channel of 2.8mm or larger before insertion.
3. Remove the device from the package and uncoil the entire device and drape in a "U" shaped configuration, holding the proximal end in one hand and the distal sheath in the opposite hand.
4. Move the finger rings back and forth to confirm that the snare loop opens and closes smoothly prior to inserting into the endoscope.
5. Ensure the snare loop is fully retracted into the catheter prior to insertion into the endoscope.

Directions for Use:

1. Insert the device into the biopsy port of the endoscope. Short strokes, 1"-1.5" (2.5cm – 3.8cm) in length, are recommended throughout device passage to avoid sheath kinking.
2. Position patient return electrode and connect to electrosurgical unit following instructions from the electrosurgical unit manufacturer.
3. Attach the respective Olympus style active cord to the snare's diathermic handle connection just prior to excision of polyps or tissue.
4. Consult the electrosurgical generator manufacturer's instructions for use for proper settings and use of the electrical surgical generator.
5. When the polyp has been endoscopically visualized, extend the distal tip of the catheter into the endoscopic field of view.
6. Gently deploy the snare loop to a fully open position.
7. Place the snare around the polyp tissue being resected using proper endoscopic technique
8. **For Rotator™ snare only:** Rotation may be accomplished by twisting the handle in the direction of the desired movement. Rotation is best achieved with the loop in the center of the lumen and not lying against the gastrointestinal wall.
9. Resect the polyp using proper surgical technique. **The electrosurgical generator should be set to the proper setting that will allow resection with controlled hemostasis.**
10. After the polyp(s) has been satisfactorily resected, the snare should be retracted back into the catheter.
11. Once the polyp(s) has been satisfactorily resected, the polyp fragment(s) should be removed and the specimen(s) prepared according to standard technique for histologic evaluation.

Product Disposal:

After use, this product may be a potential biohazard which presents a risk of cross-contamination. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Issued Date: September 2019

Warning:

An issued or revision date for these instructions is included for the user's information. In the event that two years have elapsed between this date and product use, the user should contact STERIS to determine if additional information is available.

Unless otherwise indicated, all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or are trademarks owned by STERIS Corporation.

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

US Endoscopy, a wholly owned subsidiary of STERIS Corporation.

Olympus® is a registered trademark of Olympus Corp.

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Sauf indication contraire, toutes les marques accompagnées des symboles ® ou ™ sont des marques déposées auprès de l'U.S. Patent and Trademark Office, ou sont des marques commerciales détenues par STERIS Corporation.
Sofern nicht anders angegeben, sind alle mit © oder ™ versehenen Marken beim United States Patent and Trademark Office (Patentamt der Vereinigten Staaten) registriert oder sind Eigentum der STERIS Corporation.
A menos que de outra forma indicado, todas as marcas com ® ou ™ são registadas na U.S. Patent and Trademark Office ou são marcas registadas propriedade da STERIS Corporation.
Medmindre andet er angivet er alle mærker markeret med ® eller ™ registreret hos U.S. Patent and Trademark Office eller er varemærker ejet af STERIS Corporation.
Tenzij anders is aangegeven, zijn alle aanduidingen met ® of ™ gedeponerd bij het Patent and Trademark Office in de Verenigde Staten of zijn deze aanduidingen handelsmerken van STERIS Corporation.
별도로 명시되어 있지 않는 한, ® 또는 ™으로 표시된 모든 상표는 미국 특허청(U.S. Patent and Trademark Office)에 등록되어 있거나 STERIS Corporation이 소유하고 있는 상표입니다.