# STORZ

**KARL STORZ — ENDOSKOPE**

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**INSTRUCTION MANUAL**  
IMAGE1 S H3-Z FI Camera Head TH 102

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## 1 Important note for users of KARL STORZ devices and instruments

It is recommended to check the suitability of the product for the intended procedure prior to use.

Thank you for your expression of confidence in the KARL STORZ name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high quality piece of equipment from KARL STORZ.

### 1.1 Intended use

The camera head is used together with the camera control unit and is connected to a telescope during endoscopic or microscopic interventions. It can also be used for fluorescence applications.

### 1.2 Components of the KARL STORZ FI System

The KARL STORZ FI System comprises the KARL STORZ D-LIGHT C/AF/P system including suitable KARL STORZ fluid light cables, KARL STORZ FI endoscope or FI videoscope and KARL STORZ FI camera with coordinated KARL STORZ camera control unit.

Refer to the system description for KARL STORZ Fluorescence Imaging (FI) Systems (see instruction manual 96206144DER) for the recommended configuration.

**i NOTE:** Please observe the instruction manual for the KARL STORZ camera control unit IMAGE1 S (Art. No. 96206286DER) or the instruction manual IMAGE 1 HUB™ HD (Art. No. 96206114DER).

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<b>1</b>	<b>Important note for users of KARL STORZ devices and instruments</b>	<b>III</b>
1.1	Intended use	III
1.2	Components of the KARL STORZ FI System	III
<b>2</b>	<b>Safety instructions</b>	<b>3</b>
2.1	Explanation of warnings and cautions	3
2.2	General safety information	4
<b>3</b>	<b>Images of the equipment</b>	<b>6</b>
<b>4</b>	<b>Symbols employed</b>	<b>7</b>
4.1	Symbols for operation	7
4.2	Label and packaging symbols	7
<b>5</b>	<b>Operating instructions</b>	<b>8</b>
5.1	Frequency of the High Definition camera head (50/60 Hz)	8
5.2	Connecting an endoscope and light cable	8
5.3	Adjusting the focus	9
5.4	Adjusting the image size with Zoom camera heads	9
5.5	Operating the user control buttons	10
<b>6</b>	<b>Technical description</b>	<b>11</b>
6.1	Specifications	11
6.2	Storage and operating conditions	11
6.3	Standard compliance	11
6.4	Directive compliance	11
6.5	Information on electromagnetic compatibility (EMC)	12
<b>7</b>	<b>Reprocessing</b>	<b>13</b>
7.1	Care and handling	13
7.2	Reprocessing	14
7.2.1	Reference chart	14
7.2.2	General warnings	15
7.3	Preparation for cleaning and disinfection	17
7.4	Manual wipe-down disinfection	17
7.5	Manual cleaning	18
7.6	Manual disinfection	18
7.7	Inspection and care	19
7.8	Packaging systems	19
7.9	Sterilization	19
7.10	High-Level disinfection	22
7.11	Reprocessing limit	23

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<b>8</b>	<b>Maintenance.....</b>	<b>24</b>
8.1	Maintenance and safety check.....	24
8.1.1	Maintenance.....	24
8.1.2	Safety check.....	24
8.2	Servicing and repair.....	25
8.3	Disposal.....	25
<b>9</b>	<b>Repairs and warranty.....</b>	<b>26</b>
9.1	Repair program.....	26
9.2	Limitation of liability.....	27
9.3	Warranty.....	27
<b>10</b>	<b>Compatible camera control units.....</b>	<b>28</b>
<b>11</b>	<b>Spare parts, recommended accessories.....</b>	<b>29</b>
<b>12</b>	<b>Subsidiaries.....</b>	<b>30</b>

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## 2 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

### 2.1 Explanation of warnings and cautions

Please read this manual carefully and follow its instructions precisely. The words Warning, Caution and Note convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the system. To make these words stand out more clearly, they are accompanied by a pictogram.

**WARNING:** A Warning indicates that the personal safety of the patient or physician may be at risk.

**CAUTION:** A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

**NOTE:** A Note indicates special information about operating the device or clarifies important information.

This instruction manual is intended to serve only as a supplement to the KARL STORZ video systems instruction manual, which was enclosed with your camera control unit. To ensure the safe and efficient use of this KARL STORZ product, it is essential that you read the entire video system instruction manual before connecting the camera system.

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**2.2 General safety information**

**WARNING:** Read this instruction manual thoroughly before using the equipment.

**WARNING:** Check that the equipment is functioning properly prior to each surgical intervention. Ensure that the correct video image appears on all video monitors before beginning each intervention. It is recommended that a replacement system is kept available.

**WARNING:** Perceived electric shocks may be experienced when using the camera system in conjunction with electro-surgical units.

**WARNING:** Electrostatic discharge (ESD) may be perceived as an electric shock when touching the camera system for the first time. In the event of image loss due to ESD, switch the camera control unit off and back on again to restore the image.

**WARNING:** When using high frequency surgical accessories, care must be taken to avert contact between the HF surgical device and the endoscope or endoscopic camera. Inadvertent contact may represent a thermal hazard to the surgeon.

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**WARNING:** All endoscopic devices (e.g., HF surgical probes) used in conjunction with KARL STORZ endoscopic cameras must comply with IEC/EN 60601-2-18, the standard for Medical electrical equipment – Particular requirements for the safety of endoscopic equipment. The insulation of the endoscopic devices must be designed to withstand the highest HF voltage encountered during surgery. HF-generating equipment is dangerous and cannot be allowed to make contact with the camera head or system via the accessories or endoscopes connected to the camera head or system.

**WARNING:** Risk of blindness. Never look into the open end of a light cable.

**WARNING:** Risk of burns. Contact with the open end of a light cable connected to a light source may cause burns.

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- ① Instrument coupler
- ② Focusing ring/knob
- ③ Zoom ring
- ④ User control button: for scrolling (downwards) through the menu or activating a camera function
- ⑤ User control button: for scrolling (upwards) through the menu or activating a camera function
- ⑥ User control button: for accessing the menu/ making a selection

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## 5 Operating instructions

### 5.1 Frequency of the High Definition camera head (50/60 Hz)

The IMAGE1 S H3-Z FI camera head operates on both 50 Hz and 60 Hz frequency. The frequency can be adapted to suit local conditions.

**i** **NOTE:** The method for switching frequency is described in the instruction manual IMAGE1 S – TC 200, TC 300, TC 301 (96206286D), or in the instruction manual IMAGE 1 HUB™ HD (96206114D).

### 5.2 Connecting an endoscope and light cable

The IMAGE1 camera heads have an integrated instrument coupler and optics. To connect an endoscope to the instrument coupler, turn the outer ring of the instrument coupler clockwise and insert the eyepiece of the endoscope into the coupler. When connecting the endoscope to a Pendulum camera head, the Pendulum lock must click into place before the outer ring of the grasping mechanism is turned. To prevent rotation of the endoscope, the outer ring of the instrument coupler must be turned counterclockwise until tight. Connect the light cable to the eyepiece of the endoscope by tightening the knurled screw on the light cable through a quarter turn.

**!** **CAUTION:** Heat from the light source can cause disinfectant residue on the light cable or on the light inlet of the endoscope to burn. This can have a detrimental effect on the light conduction and thus the image quality.

**i** **NOTE:** For further information regarding light cables, please refer to the manual “Fiberoptic light cables” (96206342ZS), which can be downloaded or requested at [www.karlstorz.com](http://www.karlstorz.com).

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### 5.3 Adjusting the focus

The focus can be adjusted by turning the focusing ring/knob on the camera lens.

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
### 5.4 Adjusting the image size with Zoom camera heads

The image size can be adjusted by turning the zoom ring.

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### 5.5 Operating the user control buttons

- ④ Scrolling downwards: Scroll downwards through menus, options, etc. Can also be assigned a programmable function/feature.
- ⑤ Scrolling upwards: Scroll upwards through menus, options, etc. Can also be assigned a programmable function/feature.
- ⑥ Menu/select: Access the menu. Once the menu is displayed, this button can be used to select the highlighted option (identical to Enter on a computer keyboard). Hold this button down for seven seconds to call up the “Options” menu.

For further instructions on the programming and operation of the camera head buttons, please refer to the KARL STORZ IMAGE1 S (Art. No. 96206286DER) or IMAGE 1 HUB™ HD (Art. No. 96206114DER) instruction manuals.

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**6.2 Storage and operating conditions**

Refer to the compatible camera control unit (CCU) manual for relevant storage and operating conditions. The CCU and camera head were evaluated as a system for IEC 60601-1.

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**6.3 Standard compliance**

According to IEC 60601-1; IEC 60601-2-18; UL 60601-1

- Protection of applied part against electric shock: Type CF 

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**6.4 Directive compliance**

According to the Medical Device Directive (MDD):

This medical device belongs to Class I

This medical device bears a CE mark in accordance with MDD 93/42/EEC.

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**6.5 Information on electromagnetic compatibility (EMC)**

Refer to the Electromagnetic Compatibility (EMC) information Appendix of the compatible camera control unit (CCU) manual for relevant electromagnetic compatibility (EMC) information.

The EMC warning statements, cautions, notes and emission/immunity limits specified in the EMC tables of the compatible CCU apply to the camera head. The CCU and camera head were tested as a system.

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## 7 Reprocessing

The following chapter provides information on the care and cleaning of the IMAGE1 S H3-Z FI camera head.

### 7.1 Care and handling

To maximize the life of your camera, please follow the guidelines given below for the correct handling and care of the camera head:

DO NOT PLUG THE CAMERA HEAD CABLE INTO THE PROCESSOR RECEPTACLE IF IT IS MOIST OR NOT CLEAN. Always ensure that the connector is completely clean. Remove any moisture present with a sterile cloth before connecting.

DO NOT PICK UP THE CAMERA HEAD BY ITS CABLE. Always pick up the camera head first, and then pick up and support the cable and connector. The cable contains many tiny conductors and must therefore be handled with care.

DO NOT STORE THE CAMERA SYSTEM in direct sunlight or at high temperatures.

ENSURE THAT THE CABLE IS NOT COILED TOO TIGHTLY OR FOLDED UP WHILE THE CAMERA IS IN STORAGE. Always arrange the cable in loose coils with a diameter of 6 inches (15 cm) or more.

DO NOT PULL FORCEFULLY ON THE CABLE WHEN UNCOILING. Undesirable “kinks” may result. Always allow the cable to uncoil gently and loosely.

HANDLE THE SYSTEM WITH CARE, avoiding any strong impacts. When transporting, use the original packaging to prevent damage to the device.

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● Sterilization method has been validated.

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**7.2.2 General warnings**

**WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual “Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments” and the accompanying documentation.

**WARNING:** Risk of infection: These medical devices are not sterile upon delivery. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

**WARNING:** When carrying out any work on contaminated medical devices, national regulations regarding personal safety must be observed.

**CAUTION:** When preparing and using solutions, follow the chemical manufacturer’s instructions paying close attention to proper concentration, exposure time, and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

**CAUTION:** Danger of damage to medical devices: Medical devices may be damaged by using chemicals that have not been approved by KARL STORZ.

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**CAUTION:** National laws and regulations must be observed.

**NOTE:** The manual “Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments” can be downloaded or requested by visiting [www.karlstorz.com](http://www.karlstorz.com).

**CAUTION:** Inspect the camera head cable for breaks and cracks before and after reprocessing. Camera heads with damaged cables must NOT be processed. A break or crack in the cable will allow fluids to enter the camera cable and cause damage. Return the damaged camera head to KARL STORZ for repair.

**CAUTION:** Chemical residue on the camera head plug can cause failure of the camera electronics.

**CAUTION:** The camera head must be thoroughly cleaned and disinfected/sterilized using validated infection control procedures before use and before every subsequent application. Any deviations from the recommended parameters for cleaning, disinfection and sterilization have to be validated by the user.

**CAUTION:** The use of sterilization/disinfection methods other than those recommended by KARL STORZ may have adverse effects on the materials of the camera and endoscopic adaptor. Any resulting damage will not be covered by the warranty.

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### 7.3 Preparation for cleaning and disinfection

Immerse the camera head in a suitable container containing a pH-neutral (pH value of 6.0-8.0) enzymatic cleaning solution as per the manufacturer's instructions immediately after use to prevent blood, protein and other contaminants from drying onto the equipment.

**CAUTION:** Inspect the camera head cable for breaks and cracks. Camera heads with damaged cables must not be immersed. Following wipe-down disinfection, return camera heads with damaged cables to KARL STORZ for repair.

Heavy soiling must be removed from the camera head immediately after use. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

### 7.4 Manual wipe-down disinfection

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material incompatibility must be observed.

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### 7.5 Manual cleaning

The camera head and camera head cable must be completely immersed, preferably in a pH-neutral enzymatic cleaning solution or microbiologically pure/completely demineralized water. It must be ensured that the surfaces are covered and that no air bubbles are present. After the required exposure time, clean the instrument mechanically with a soft brush (Art. No. 27652) or sponge. Finally, it must be rinsed with microbiologically pure, sterile or completely demineralized water to ensure neutralization. Thorough rinsing of the camera head is necessary for removing any soiling and cleaning agent residues which could interfere with disinfection.

### 7.6 Manual disinfection

The camera head and the camera head cable must be completely immersed in a disinfectant solution. It must be ensured that the surfaces are covered and that no air bubbles are present. Following exposure, rinse the camera head three times with microbiologically pure water, sterile water, or completely demineralized water. Each rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Thorough rinsing of the camera head is necessary for removing any residual disinfectant that could interfere with sterilization.

All the surfaces are then dried completely with sterile compressed air. The cleaning gun with accessories (Art. No. 27660) is ideal for this purpose. Wipe the camera head connector with a cotton tip applicator moistened with 70% isopropyl alcohol to remove any residual disinfectant.

Do not allow exposed glass windows to air dry. 70% isopropyl alcohol may be applied to glass surfaces with a cotton tip applicator to prevent streaks and spots. Dry the surfaces thoroughly with a cotton tip applicator after applying the alcohol.

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### 7.7 Inspection and care

The cleaned camera head and cleaned camera head cable must be visually inspected for cleanliness, completeness, damage and dryness.

- If residues or contamination are still present, the camera head and the camera head cable must be manually cleaned and subjected to a full cleaning procedure once more.
- Camera heads with damaged cables must be withdrawn from use.

**CAUTION:** Inspect the camera head cable for breaks and cracks. Camera heads with damaged cables must not be sterilized or disinfected. Following wipe-down disinfection, return camera heads with damaged cables to KARL STORZ for repair.

- Afterwards, a functional check must be carried out.

Before sterilization and/or disinfection, coil the camera head cable into loops with a diameter of at least 15 cm (6"). In doing so, do not kink or twist the cable.

### 7.8 Packaging systems

Only standardized and approved packaging materials or systems may be used (EN 868 Parts 2 – 10, EN ISO 11607 Parts 1 + 2, DIN 58953). National regulations must be observed.

### 7.9 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments". The procedure must be selected taking the respective applicable national requirements into account and in consultation with the device manufacturers.

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The following sterilization methods have been validated and approved by KARL STORZ for the IMAGE1 S H3-Z FI camera head:

**Hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) sterilization – ASP STERRAD®**

**WARNING:** Remember that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of dimensions and material.

**CAUTION:** For more detailed information, please consult the user handbook of the respective device.

**NOTE:** To ensure that the medical device concerned can be sterilized using the different STERRAD® units, refer to the “STERRAD® Sterility Guide”.

**CAUTION:** STERRAD® sterilization may cause cosmetic changes to the camera head usually do not impact the camera's functionality.

The following STERRAD® sterilization methods have been validated and approved by KARL STORZ for the IMAGE1 S H3-Z FI camera head:

- STERRAD® 100S „Short“ cycle
- STERRAD® NX „Standard“ cycle
- STERRAD® 100NX „Standard“ cycle
- „DUO“ cycle

**Hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) sterilization – STERIS® AMSCO®**

More detailed information regarding the cycle used for various device generations can be obtained from the manufacturer STERIS®.

**CAUTION:** V-PRO™ sterilization may cause cosmetic changes to the camera head usually do not impact the camera's functionality.

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The following STERIS® AMSCO® sterilization methods have been validated and approved by KARL STORZ for the IMAGE1 S H3-Z FI camera head:

- V-PRO™ 1  
Equivalent to “lumen” cycle
- V-PRO™ 1 Plus  
“Lumen” cycle  
“Non-Lumen” cycle
- V-PRO™ maX  
“Lumen” cycle  
“Non-Lumen” cycle  
“Flexible” cycle

**Ethylene oxide sterilization (EO)**

The ethylene oxide procedure is validated with 100% ethylene oxide at 55°C and a hold time of 30-45 minutes.

**i NOTE:** Place the camera head in a sterilization tray. Do not place any instruments on top of the camera.

**! CAUTION:** For gas sterilization with ethylene oxide, the airing times stipulated for the materials by the device manufacturer must be observed due to gas absorption by the materials. The airing times depend on the processing procedure of the EO sterilizer (concentration, process control).

**Chemical low-temperature sterilization with peracetic acid – STERIS® System 1 (except U. S.)**

More detailed information regarding sterilization parameters can be obtained from the manufacturer STERIS® .

**Chemical low-temperature sterilization with peracetic acid - STERIS® System 1E (only valid in U. S.)**

More detailed information regarding sterilization parameters can be obtained from the manufacturer STERIS® .

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**7.10 High-Level Disinfection**

**WARNING:** High-level disinfection should ONLY be used for instruments which come into contact with intact mucous membranes.

**CAUTION:** Any deviations from the recommended parameters for disinfection must be validated by the user.

**CAUTION:** Before disinfection, the instruments must be thoroughly cleaned, rinsed, and dried.

**CAUTION:** When preparing and using solutions, follow the chemical manufacturer's instructions paying close attention to proper concentration, exposure time, and service life.

The following high-level disinfection methods have been validated and approved by KARL STORZ for the IMAGE1 S H3-Z FI camera head:

**CIDEX (2.4% glutaraldehyde solution)**

In order to prevent potential damage, do not place camera head in a container with other instruments. Completely immerse the camera head in the undiluted 2.4% glutaraldehyde solution for a minimum of 45 minutes at 25°C (77°F). It must be ensured that the surfaces are covered and that no air bubbles are present.

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To finish, rinse the camera head three times with sterile water. Each rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Use fresh sterile water for each rinse. Thorough rinsing of the camera head and the card-edge connector with sterile water is very important for removing any residual disinfectant.

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Dry the camera head and the card-edge connector with a lint-free, sterile cloth. Ensure that the card-edge connector is completely clean and free of moisture before connecting it to the camera processor. Wipe the camera head connector with isopropyl alcohol to remove any residual disinfectant.

**Resert XL HLD (2.0% hydrogen peroxide solution)**

In order to prevent potential damage, do not place camera head in a container with other instruments. Completely immerse the camera head in the undiluted 2.0% hydrogen peroxide solution for a minimum of 8 minutes at 20°C (68°F). It must be ensured that the surfaces are covered and that no air bubbles are present.

**CAUTION:** Do not leave the camera head in the solution for any longer than 60 minutes.

To finish, rinse the camera head three times with sterile water. Each rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Use fresh sterile water for each rinse. Thorough rinsing of the camera head with sterile water is very important for removing any residual disinfectant.

Dry the camera head and the card-edge connector with a lint-free, sterile cloth. Ensure that the card-edge connector is completely clean and free of moisture before connecting it to the camera processor. Wipe the camera head connector with isopropyl alcohol to remove any residual disinfectant.

**7.11 Reprocessing limit**

The end of the product's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

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### 8.1 Maintenance and safety check

#### 8.1.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be requested from your local representative or from the manufacturer.

#### 8.1.2 Safety check

Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, the following checks must be carried out on this device at least once a year by persons who, due to their training, knowledge and practical experience, are able to properly carry out such safety checks and who are not subject to any instructions as far as such checking activities are concerned.

##### Visual inspection:

- Check the device and accessories for any mechanical damage which may impair functionality.
- Check that the inscriptions relevant to safety are legible.

##### Electric measurements:

- Protective grounding measured according to IEC 62353 (current edition), and checked for compliance with the limits specified in the standard.
- Touch Current, Earth Leakage Current, and Patient Leakage currents measured according to IEC 62353 (current edition), and checked for compliance with the limits specified in the standard.

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**Test for proper functioning:**

- A functional check must be carried out in accordance with the instruction manual (see Chapter 7.3.12 of the IMAGE1 S CCU instruction manual no. 96206286 DER).
- The safety check must be entered in the unit log-book and the results of the check must be documented.
- If the device does not function reliably and/or safely, it must be repaired.

**8.2 Servicing and repair**

Defective items of equipment must be serviced and repaired exclusively by KARL STORZ or by persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

**CAUTION:** Do not attempt to open the camera head or lens. Opening these sealed assemblies affects the camera's soakability and voids all warranties. Do not remove the covers on the camera processor. Warranty is voided unless removed by authorized personnel.

**8.3 Disposal**

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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**9 Repairs and warranty****9.1 Repair program**

**CAUTION:** Opening the camera head will result in damage and void the warranty.

To bridge the repair period, you will, as a rule, receive a unit on loan which must then be returned to KARL STORZ as soon as you receive the repaired unit.

In the case of repairs, please contact your local KARL STORZ subsidiary or authorized dealer.

In Germany you can refer repairs direct to  
KARL STORZ GmbH & Co. KG  
Abt. Reparaturservice  
Take-off Gewerbepark 83  
78579 Neuhausen, Germany

**Important information**

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among personnel).

We reserve the right to return contaminated instruments/devices to the sender.

Repairs, modifications, or expansions which are not performed by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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### 9.2 Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the unit is connected and operated meet the applicable laws and standards; and
- the unit has been used in accordance with its operating instructions at all times.

### 9.3 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see “Subsidiaries” section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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**i NOTE:** The S technologies are available with the IMAGE1 S.

\* IMAGE1 camera control units are available with various configurations and image outputs. These camera control units are compatible with IMAGE1 S camera heads.

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**11 Spare parts, recommended  
accessories**

Article	Order no.
Adaptor, autoclavable, allows telescopes to be changed under sterile conditions	533 TVA
Integral rotating beamsplitter, for simultaneous viewing through the telescope and on the monitor	9530 BD
Sterile drapes* for KARL STORZ cameras	
Extension cable 5.2 m	22 2200 72
Dust cap for camera heads	6349190
Plastic container for sterilization and storage of IMAGE1 S H3-Z FI (TH102) camera head	39301 Z3TS
Instruction manual	96206520DER

\* For safe and comfortable working with KARL STORZ cameras, we recommend the use of sterile drapes.

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KARL STORZ Endoscopy Canada Ltd.  
7171 Millcreek Drive, Mississauga  
Ontario, L5N 3R3, Canada  
Phone: +1 905 816-4500, Fax: +1 905 858-4599  
Toll free phone: 1-800-268-4880 (Canada only)  
Toll free fax: 1-800-482-4198 (Canada only)  
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.  
2151 East Grand Avenue  
El Segundo, CA 90245-5017, USA  
Phone: +1 424 218-8100, Fax: +1 424 218-8526  
Toll free phone: 800 421-0837 (USA only)  
Toll free fax: 800 321-1304 (USA only)  
E-Mail: info@kseas.com

KARL STORZ Veterinary Endoscopy-America, Inc.  
175 Cremona Drive  
Goleta, CA 93117, USA  
Phone: +1 805 968-7776, Fax: +1 805 685-2588  
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.  
815 N. W. 57th Avenue, Suite 480  
Miami, FL 33126-2042, USA  
Phone: +1 305 262-8980, Fax: +1 305 262-8986  
E-Mail: info@ksela.com

KARL STORZ Endoscopia  
Miramar Trade Center,  
Edificio Jerusalem, Oficina 108,  
La Habana, Cuba  
Phone: +53 72041097, Fax: +53 72041098  
E-Mail: karlstorz@enet.cu

KARL STORZ Endoscopia México S.A de C.V.  
Lago Constanza No 326,  
Col. Granada Del. Miguel Hidalgo  
C.P. 11520 México D.F.  
Phone: +52 (55) 1101 1520  
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.  
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,  
CEP-04534-011 São Paulo, Brasil  
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680  
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.  
Zufriategui 627 6º Piso, B1638 CAA - Vicente Lopez  
Provincia de Buenos Aires, Argentina  
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773  
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS  
Østensjøveien 15B  
0661 Oslo, Norway  
Phone: +47 6380 5600, Fax: +47 6380 5601  
post@karlstorz.no

KARL STORZ Endoskop Sverige AB  
Storsåtragränd 14  
12739 Skärholmen, Sweden  
14108 Kungens Kurva, Sweden  
Phone: +46 8 505 648 00  
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY  
Valimotie 1 B, 4. kerros  
00380 Helsinki, Finland  
Phone: +358 (0)96824774, Fax: +358 (0)96824775  
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ GmbH & Co. KG  
Representation Office  
Kestucio g. 59  
08124 Vilnius, Lithuania  
Phone: +370 68567000  
E-Mail: ausra.kukeniene@karlstorz.com

KARL STORZ Endoskopi Danmark A/S  
Skovlytoften 33  
2840 Holte, Denmark  
Phone: +45 45162600, Fax: +45 45162609  
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.  
415 Perth Avenue, Slough  
Berkshire, SL1 4TU, United Kingdom  
Phone: +44 1753 503500, Fax: +44 1753 578124  
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B. V.  
Displayweg 2  
3821 BT Amersfoort, Netherlands  
Phone: +31 (0)33 4545890  
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.  
Phone: +32 473 810 451  
E-Mail: info@stopleer.be

KARL STORZ Endoscopie France S. A. S.  
12, rue Georges Guynemer, Quartier de l'Europe  
78280 Guyancourt, France  
Phone: +33 1 30484200, Fax: +33 1 30484201  
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH  
Landstraßer Hauptstr. 148/1/G1  
1030 Wien, Austria  
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479  
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.  
Parque Empresarial San Fernando  
Edificio Munich – Planta Baja  
28830 Madrid, Spain  
Phone: +34 91 6771051, Fax: +34 91 6772981  
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.  
Via dell'Artigianato, 3  
37135 Verona, Italy  
Phone: +39 045 8222000, Fax: +39 045 8222001  
E-Mail: info-ita@karlstorz.com

KARL STORZ Adria Eos d.o.o.  
Capraška 6  
10000 Zagreb, Croatia  
Phone: +385 1 6406 070, Fax: +385 1 6406 077  
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.  
Verovškova c. 60A  
1000 Ljubljana, Slovenia  
Phone: +386 1 620 5880, Fax: +386 1 620 5882  
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.  
ul. Bojkowska 47  
44-100 Gilwice, Poland  
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07  
E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy India Private Limited  
11th Floor, Dr. Gopal Das Bhawan  
28, Barakhamba Road  
New Delhi 110001, India  
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010  
E-Mail: corporate@karlstorz.in

KARL STORZ GmbH & Co. KG  
Resident Representative Office  
14th Floor, Kumho Asiana Plaza  
39 Le Duan, District 1  
Ho Chi Minh City, Vietnam  
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039  
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.  
15 Orion Road Lane Cove NSW 2066  
P O Box 50 Lane Cove NSW 1595, Australia  
Phone: +61 (0)2 9490 6700  
Toll free: 1800 996 562 (Australia only)  
Fax: +61 (0)2 9420 0695  
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong  
Room 2503-05, 25F AXA Tower, Landmark East,  
No. 100 How Ming Street, Kwun Tong, Kowloon,  
Hong Kong, People's Republic of China  
Phone: +852 28 65 2411, Fax: +852 28 65 4114  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.  
Beijing Branch,  
Room 1805-1807, Building B, 18F Beijing IFC  
No. 8, Jianguomenwai Street, Chaoyang District,  
Beijing, 100022, People's Republic of China  
Phone: +86 10 5638188, Fax: +86 10 5638199  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Shanghai Branch  
36F, Tower 1 Grand Gateway  
No. 1 Hongqiao Road, Xuhui District,  
Shanghai, 200030, People's Republic of China  
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Chengdu Branch  
Room 805, 8F Jin Jiang International Building  
No. 1 West Linjiang Road, Wuhou District,  
Chengdu, 610044, People's Republic of China  
Phone: +86 28 86587977, Fax: +86 28 86587975  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Shenyang Branch,  
Room 2001-2005, 20F N-MEDIA International Center,  
No. 167 Youth Avenue, Shenhe District,  
Shenyang, 110014, People's Republic of China  
Phone: +86 24 23181118, Fax: +86 24 23181119  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Guangzhou Branch  
Room 02B, 35F Teem Tower,  
No. 208 Tianhe Road, Tianhe District,  
Guangzhou, 510620, People's Republic of China  
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.  
No. 8 Commonwealth Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd  
No. 8 Commonwealth Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.  
I-Park Bundang1, 102dong, 8F  
239 Jeongjail-Ro, Bundang-Gu  
Seongnam, 463-859, Korea  
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299  
E-Mail: stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.  
12F, No. 192, Sec. 2, Chung Hsin Rd.,  
Sindian District, New Taipei City, Taiwan  
Phone: +886 933 014 160, Fax: +886 2 8672 6399  
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.  
Stage Bldg. 8F, 2-7-2 Fujimi  
Chiyoda-ku, Tokyo 102-0071, Japan  
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633  
E-Mail: info@karlstorz.co.jp

www.karlstorz.com

**STORZ**  
KARL STORZ—ENDOSKOPE

**KARL STORZ GmbH & Co. KG**

Mittelstrasse 8  
78532 Tuttlingen

Postfach 230  
78503 Tuttlingen  
Germany

Telefon: +49 (0)7461 708-0  
Telefax: +49 (0)7461 708-105  
E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)  
Web: [www.karlstorz.com](http://www.karlstorz.com)