

300 Series

Infinium Luxor 300 Mobile Lighting Platform

A sleek high adjustable intensity LED surgical lighting system with touch screen and wireless remote features. This well built line is designed to minimize contamination from blood, chemicals and air circulation and is easy to clean.



Illumination Intensity
3 Step adjustment (Max. 130,000 lx)



Light Field Size
3 Step adjustment (Max. 9.84in / Min. 5.9in)



Color Temperature
3,800K / 4,300K / 4,800K



Color Rendering Index
Ra 96 / R9 95



3 Different Modes
Surgical mode (130,000 lx)
Examine mode (60,000 lx)
Dental mode (30,000 lx)



Optional Configurations

