

300 Series

300 Series

Advanced Surgical Lighting

A sleek high adjustable intensity LED surgical lighting system with touch screen and wireless remote features. This well built line is designed to minimize contamination from blood, chemicals and air circulation and is easy to clean.



Illumination Intensity

3 Step adjustment (Max. 130,000 lx)



Light Field Size

3 Step adjustment (Max. 9.84in / Min. 5.9in)



Color Temperature

3,800K / 4,300K / 4,800K



Color Rendering Index

Ra 96 / R9 95



3 Different Modes

Surgical mode (130,000 lx)

Examine mode (60,000 lx)

Dental mode (30,000 lx)




Optional Configurations



Bright, Reliable, & Intuitive

TECHNICAL SPECIFICATIONS

Class	S300
Design (Light Head)	
Installation	Single Ceiling Dual Ceiling Mobile
Head Size (in)	17.1 in
No. of LED (ea)	42
Intensity (lx)	130,000
Intensity Control	3 Steps (40,000~130,000)
C.R.I (Ra/R9)	96 / 95
Color Temp. (K)	3,800 / 4,300 / 4,800
Light Field Size (in)	3 Steps (5.9~9.84, Circle)
Detachable Centre Grip	Yes (2X)
Input Power	100~240Vac, 50/60Hz
LED Life(Hrs)	60,000
Option	Wall Controller (Tablet PC) Backup Battery

Battery Backup


- Over & Under Voltage Detection
- Discharging Over Current Detect
- Short Circuit Detection
- Cell Balancing, PTC
- 300 Cycle Rating
- 3 Hour Backup Charge Duration

Specifications	
Rating	20~28Vdc, 4A (Lithium-ion Battery)
Dimension	9.6 x 6.4 x 4.7in
Weight	8.38 lbs

300 Series Surgical Lights

Premium surgical, operating room, and procedure lights at a price to fit your budget. The 300 series is long lasting, feature rich, easy to use and backed by a 3yr warranty

Modern Touch Control Systems



OR Light (1 or 2) Control

- Light 1/ Light 2 Select
- Light ON/OFF
- Field Size, Intensity, Color Temp

Sync Mode

- OR Light (1 and 2) Sync

Dental Mode

- Dental Light ON/OFF
- Intensity Control

Max. Mode OR Light

- Light Field Size 3 Level and Intensity 3 Level Setting

Examination Mode

- Exam. Light ON/OFF
- Intensity, Color Temp. Control


Light On/Off

Hand Touch Controls

Intensity Control

(1 → 2 → 3 → 1, Rotate)

Field Adjustment



INFINIUM®

Infinium Medical
12151 62nd St North #5
Largo, FL 33773 • USA

(1) 866-918-8434
infiniummedical.com
sale@infiniummedical.com