




Clip hemostático desechable



disposable

Back

Disposable HemoClip



Item No.	Tube Diameter & Working Length	Working Channel Diameter	Opening Width	Scope Of Use
AF-D2418JZR-9	Ø2.4x1800mm	≥2.8mm	9mm	Gastroscopy
AF-D2423JZR-9	Ø2.4x2300mm	≥2.8mm	9mm	Colonoscopy
AF-D2418JZR-11	Ø2.4x1800mm	≥2.8mm	11mm	Gastroscopy

Rotate 360° in different directions. Large opening design, Meet clinical treatment needs. Large opening design, Meet clinical treatment needs. Short tail, Easy to release, More safety. Opening Width: 9mm / 11mm / 14mm / 17mm