



Aguja de inyección desechable



Disposable Injection Needle



Item No.	Tube Diameter & Working Length	Working Channel Diameter	Needle Length	Needle Diameter	Scope Of Use
AF-D1816PN	Ø1.8x1600mm	≥2.0mm	4/5/6mm	23/25G	Gastroscopy
AF-D2416PN	Ø2.4x1600mm	≥2.8mm	4/5/6mm	21/23/25G	Gastroscopy
AF-D2421PN	Ø2.4x2100mm	≥2.8mm	4/5/6mm	21/23/25G	Colonoscopy
AF-D2423PN	Ø2.4x2300mm	≥2.8mm	4/5/6mm	21/23/25G	Colonoscopy

For cutting small polyps and coagulation of small vessels.