





Dilation Balloon Catheter



Usage/ For the dilation for digestive and respiratory tract stenosis.



Characteristic

- High pressure-resistant materials ensure safe dilation.
- Cylindrical balloon is not easy to deform, and after dilation, it is convenient to observe.
- Various balloon size for choose, meets different clinical demand.

Dilation Balloon Catheter

Dilation balloon for ERCP

Model	Balloon O.D. (mm)	Balloon Length (mm)	Catheter Length (mm)	Working Pressure (atm)	Remark
VDK-BD-06-40-200	6	40	2000	5	Suitable for Endoscopic channel $\geq 2.8\text{mm}$; Guide wire diameter: 0.035"
VDK-BD-08-40-200	8	40	2000	5	
VDK-BD-10-40-200	10	40	2000	5	
VDK-BD-12-40-200	12	40	2000	5	
VDK-BD-14-40-200	14	40	2000	5	
VDK-BD-16-40-200	16	40	2000	5	

Dilation balloon for PTCS (percutaneous transhepatic cholechofibroscopy)

Model	Balloon O.D. (mm)	Balloon Length (mm)	Catheter Length (mm)	Working Pressure (atm)	Remark
VDK-BD-06-40-90	6	40	900	5	Suitable for Endoscopic channel $\geq 2.0\text{mm}$; Guide wire diameter: 0.035"
VDK-BD-08-40-90	8	40	900	5	
VDK-BD-10-40-90	10	40	900	5	Endoscopic channel $\geq 2.8\text{mm}$

Dilation balloon for Cardia

Model	Balloon O.D. (mm)	Balloon Length (mm)	Catheter Length (mm)	Working Pressure (atm)	Remark
VDK-BD-25-80-90	25	80	900	2	Not suitable for Endoscope use
VDK-BD-30-80-90	30	80	900	2	
VDK-BD-35-80-90	35	80	900	2	Guide wire with diameter 0.035" can pass through.
VDK-BD-40-80-90	40	80	900	2	

Dilation balloon for Esophagus, Pylorus and Intestinal tract

Model	Balloon O.D. (mm)	Balloon Length (mm)	Catheter Length (mm)	Working Pressure (atm)	Remark
VDK-BD-10-40-180	10	40	1800	5	Suitable for Endoscopic channel $\geq 2.0\text{mm}$; Guide wire diameter: 0.035"
VDK-BD-12-60-180	12	60	1800	5	
VDK-BD-14-60-180	14	60	1800	5	
VDK-BD-16-80-180	16	80	1800	5	
VDK-BD-18-60-180	18	60	1800	4	
VDK-BD-18-80-180	18	80	1800	4	Endoscopic channel $\geq$ 3.2mm; Not suitable for Guide wire.
VDK-BD-20-60-180	20	60	1800	4	

Other sizes

Model	Balloon O.D. (mm)	Balloon Length (mm)	Catheter Length (mm)	Working Pressure (atm)	Remark
VDK-BD-06-30-110/90	6	30	1100/900	5	Suitable for Endoscopic channel $\geq 2.0\text{mm}$; Guide wire diameter: 0.032"
VDK-BD-08-30-110/90	8	30	1100/900	5	
VDK-BD-10-40-110/90	10	40	1100/900	5	
VDK-BD-12-40-110/90	12	40	1100/900	5	Suitable for Endoscopic channel $\geq 2.8\text{mm}$; Guide wire diameter: 0.035"
VDK-BD-14-60-110/90	14	60	1100/900	5	
VDK-BD-16-60-110/90	16	60	1100/900	5	
VDK-BD-18-60-110/90	18	60	1100/900	4	