

ZKSK® Luxury PEG Kit



SPECIFICATIONS

- | | |
|--|---|
| 01. 1 Traction Removable PULL PEG Tube | 09. 1 Disposable Scalpel with #11 Blade |
| 02. 1 Feeding Adapter | 10. 1 Introducer Needle |
| 03. 1 Feeding Adapter with Luer lock | 11. 1 Tubing Clamp |
| 04. 1 Retention Ring | 12. 4 10*10 Gauze Sponges |
| 05. 1 19Ga Needle | 13. 1 Fenestrated Drape |
| 06. 1 25Ga Needle | 14. 1 Direction for Use |
| 07. 1 Placement Wire (0.6mm*260cm) | 15. 1 Scissors & 1 Tweezers |
| 08. 1 10ml Syringe | 16. 1 Polypectomy Snare |

REFERENCE NUMBER

GT12F-L、GT14F-L、GT16F-L、GT18F-L、**GT20F-L**、GT22F-L、GT24F-L



Beijing ZKSK Technology Co., Ltd

Building 9, 6 & No.6 Yuan Hengye North 7th street, Yongle Economic Development Zone,
Tongzhou District, Beijing 101105, China

Tel: 86 10 57272222

E-mail: len@bjzsk.com


www.bjzsk.com

EC REP Prolinx GmbH

Brehmstr. 56, 40239 Duesseldorf, Germany




Percutaneous Endoscopic Gastrostomy (PEG) Kit & Replacement Tube



Beijing ZKSK Technology Co., Ltd



ZKSK® Silicone Gastrostomy Tube



INDICATIONS

The Gastrostomy Feeding Tube allow for delivery of enteral nutrition and medication directly into the stomach and/or gastric decompression. Mainly suitable for Gastrostomy patients.

ADVANTAGES

- Minimizing trauma during surgery.
- Made of 100% medical grade silicone, the tube is soft & clear.
- X-ray opaque line through the whole tube.
- The balloon is glued to the main tube both inside and outside, it is elastic and flexible.
- Fully equipped, easily operated.
- Good biocompatibility.
- Y Type locking joint, no leakage.
- The size from 12Fr to 24Fr, color code for distinguish different size.

SPECIFICATIONS

Type	Reference Number	Size(Fr)	Color Code	Length (mm)	Marked Position (cm)
Tubes	GT12F	12Fr/5ml	White	240	2/4/6
	GT14F	14Fr/5ml	Green		
	GT16F	16Fr/20ml	Orange		
	GT18F	18Fr/20ml	Red		
	GT20F	20Fr/20ml	Yellow		
	GT22F	22Fr/20ml	Purple		
	GT24F	24Fr/20ml	Dark Blue		

ZKSK® Standard PEG KIT



INDICATIONS

It is used for oral feeding difficulties, but gastrointestinal function is normal, the need for long-term nutritional support, including Neurotic coma, esophageal cancer, anorexia, Alzheimer's disease, etc. Mainly suitable for Gastrostomy patients.

SPECIFICATIONS

- 01. 1 Traction Removable PULL PEG Tube
- 02. 1 Feeding Adapter
- 03. 1 Feeding Adapter with Luer lock
- 04. 1 Retention Ring
- 05. 1 19Ga Needle
- 06. 1 25Ga Needle
- 07. 1 Placement Wire (0.6mm*260cm)
- 08. 1 10ml Syringe
- 09. 1 Disposable Scalpel with #11 Blade
- 10. 1 Introducer Needle
- 11. 1 Tubing Clamp
- 12. 4 10*10 Gauze Sponges
- 13. 1 Fenestrated Drape
- 14. 1 Direction for Use



REFERENCE NUMBER

GT12F-K、GT14F-K、GT16F-K、GT18F-K、GT20F-K、GT22F-K、GT24F-K